# Laser Doppler Velocimetry

# LDA

Grab your Flow!

Intelligent Laser  
Applications GmbH

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The Laser Doppler Velocimetry (LDV) is a non-intrusive method for measuring flow velocities. LDV is known for high accuracy (up to 0.1%), a wide measuring range and for its high spatial and time resolution. The ILA is able to offer a wide range of LDV systems, starting with a reasonable priced basic system up to a 3D- LDV system, due to the modularity of the components.

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A pattern of interference fringes is produced at the crossing point of two laser beams. Tiny particles (1-10  $\mu\text{m}$ ) are backscattering the laser light with a frequency proportional to the flow velocity when passing the measuring point. The backscattered light is converted with a photodetector to an electrical signal. By using FFT technics the Doppler frequency is determined. The flow velocity is calculated with the known distance between the interference fringes and the measured frequency.

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
The measuring system consists of laser, optic, controller, FFT modul and software. A HeNe laser, which is directly coupled with the optic (fp40 or fp50), is used for the basic system as light source. The advantage of this setup compared to fiberoptical solutions is the high efficiency of the laser beam delivering and the low costs. The space-saving fiberoptical systems can be equipped with Ar<sup>+</sup> and Nd:YAG lasers. The fiberoptical fp63 probe can be easily converted into a 2D-system by adding additional transmitting fibers. The coupling of the laser light into the fibers, the wavelength separation, and the frequency shifting is done within a fiber box.

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
The evaluation of the doppler bursts by FFT methods is done with the well established PC based spectral analysis modules. Datarates of up to 1.8 kHz are achievable. The software contains a traversing module for positioning traversing units with up to three axes. An optional beam calculation module allows the determination of the measuring points due to curved windows.

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
Special assignments require often specific solutions. The ILA offers next to the standard systems also the development of custom specific solutions e.g. optics for high accuracy volume flow rate measurements. Automated wind tunnels with integrated LDV-system as a standard for calibration of other flow sensors are available.

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
150 mW Nd:YAG-Laser integrated in controller unit, also available with Ar<sup>+</sup>-laser, one transmitting fiber (monomode, polarisation preserving), Braggcell and photomultiplier integrated in probe fp80, focal length  $f = 310-500$  mm, velocity range  $v < 100$  m/s

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Fiberprobe, focal length  $f = 160-800$ mm, 2 transmitting fibers, 1 receiving fiber, the 1D-probe can be upgraded to a 2D-System, watertight, fiberbox with Braggcell and dispersion prism for wavelength separation, including high precision input manipulators, photomultiplier integrated in Controller

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Designed for underwater applications, focal length in water  $f = 600$ mm, tube for easy measuring volume adjustment

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The ILA is also producing PIV systems for measuring the flow velocity of a whole flow field at a certain time. Our know how in fluid mechanics is available as a service to help to solve your problem. We would be glad to make a demonstration in your application. We have the intelligent solution.

**Intelligent Laser Applications GmbH**

Karl-Heinz-Beckurts-Straße 13  
52428 Jülich  
Telefon: 02461. 690 430  
Telefax: 02461. 690 439  
e-mail: info@ila.de  
www.ila.de