
# PIV

## Particle Image Velocimetry

Grab your Flow!

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The Particle Image Velocimetry (PIV) is a non-intrusive laser optical method for a fast acquisition of whole flow fields. The flow vector can be determined in a measuring plane with two (2D-PIV) or all three components (stereoscopic 3D-PIV). A 3D data set of the fluid velocities can be generated by moving the measuring plane through the interesting volume. PIV is already used in all areas of fluid mechanics like aerodynamic and space technology, chemical industry, mechanical and hydraulic engineering and car industry.

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
A laser is generating a thin light sheet inside the flow. With the pulsed laser a cross correlation camera is acquiring two consecutive images of particles transported with the flow. The particles have a displacement proportional to the flow velocity. With the known pulse distance and the calibrated scaling factor of the camera the flow velocity can be calculated. The range of flow velocities starts at several mm/s and can go up to transsonic flows by adjusting the laser pulse distance.

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All the components can be mounted in a trolley. The result is a mobile, compact and powerfull measuring system, which can provide measurements in several locations in your lab. The PIV system can be extended at any time with new components like e.g. cameras for stereoscopic PIV and Laser Induced Fluorescence (LIF), backplane tilt adapters and traversing units for volume scanning, due to the modular concept of the whole system.

Mobiles PIV System / mobile PIV system

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
The PIV system consists of the light sheet generation, the synchronisation and the image acquisition and evaluation. The light sheet is typically generated by a pulsed Nd:YAG Laser. The laser power is between 25 and 350 mJ/pulse. A light arm is delivering the very powerfull laser beam to the light sheet optic. Very accurate mirrors with a special coating inside the light arm are transmitting the laser light nearly without any losses. The light sheet optic generates the light sheet out of the laser beam, where the thickness and the divergence angle is adjustable infinitely variable. The synchronisation between laser and camera is controlled by the PIV Synchronizer, which can be triggered by external events. The image acquisition and evaluation is realized by a special cross correlation camera and a very fast PC. The image data is transferred with fiber optics and is evaluated, archived and displayed by a special PIV software (VidPIV).

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A special importance has the PIV software due to the amount of acquired image data and their complexity. VidPIV ensures an easy image acquisition, puts the images in its unique tree structure (archive) and has a powerfull evaluation and postprocessing of the measuring data. VidPIV offers diverse Import and Export functions, efficient Mappings, nearly any Annotations, an own Software Development Kit and supports nearly any hardware. The numerous features of VidPIV are described in an own software brochure.

Stereoskopisches 3D-PIV System / stereoscopic 3D-PIV system

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
External PIV Synchronizer with LC display and front panel menu buttons, shortest pulse distance 125 ns, outputs for laser and camera control, external trigger input for phase correlated image acquisition

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Cooled CCD chip with 1280x1024 pixels, 200 ns shortest frame straddling time, 12 bit grey values, very high quantum efficiency, 4 measurements/s with full resolution or up to 16 with binning, image data transmitted by fiber optics

▼ Mobiles PIV System / Mobile PIV system

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Laser: Nd:YAG from 25 up to 120 mJ/pulse Light arm: working range 1600 mm, seven elements, carbon fiber design, coated mirrors Light sheet optic: variable light sheet angle btw. 10-20°, parallel 50 mm or 100 mm, endoscopic light sheet available Correlation camera: Peltier cooled, 1280x1024 pixels, minimum frame straddling time 200 ns, 12 bit grey values, very high quantum efficiency Software: VidPIV NT/2000, Stereo, LIF and Multiprocessor support

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Some PIV applications need additional support e.g. endoscopic light sheet optics, watertight casings, special light delivery systems or customized software solutions which the ILA offers. The ILA can put together a customized PIV system – even with the use of existing components. Beside PIV systems, ILA also offers Laser Doppler Anemometers for precise pointwise measurements. Our fluid mechanic know how is also available as service for solving your problems. Therefore the ILA is the ideal partner in the range of non intrusive laseroptical fluid measurements.

Intelligent Laser  
Applications GmbH

Karl-Heinz-Beckurts-Straße 13  
52428 Jülich  
Telefon: 02461. 690 430  
Telefax: 02461. 690 439  
e-mail: [info@ila.de](mailto:info@ila.de)  
[www.ila.de](http://www.ila.de)