


KL 1500 LED plus

LED light source with multiple LED light engine



Application:

illumination via fiber optics for industrial and life science applications of stereo microscopy

Main Features:

- basic features like KL 1500 LED
- remote controllable
 - PC via USB
 - foot switch
- display
 - brightness 0...100%
 - status
 - error codes
- electronic shutter function
 - flash (freezing of movements)
 - time lapse



Highlighted Technical Data:

- luminous flux: 380 lm out of an 8 mm diameter fiber bundle
- maximum power consumption: 36 W



SCHOTT
glass made of ideas